
Your TRUETAPE® Crosstapes

TRUETAPE® Crosstapes are available in 3 different sizes. They can be used almost all over the body. You can apply them directly onto pain / trigger / and acupuncture points*.

TRUETAPE® Crosstapes are waterproof, skin friendly and completely drug free. The wearing time varies between 2 and 4 days, depending on the area of application. For the best taping experience it is recommended to renew the tapes on a regular basis.

Materials: 100% viscose with hypoallergenic acrylic adhesive
In case of allergies to any of the ingredients please do not use TRUETAPE®.

How it works:

1. Locate your pain or trigger point*
2. Choose appropriate Crosstape size
3. Apply Crosstape to the skin and rub gently

The application*

Apply your TRUETAPE® Crosstapes directly to pain points* on muscles or joints. Always make sure to apply them only to healthy and clean skin.

Choose ideal tape size

All Crosstapes can be universally used almost on the entire body. However, some tape sizes are even better suited for certain body parts than others:

Size S: Ideal for kids or for the head, throat, feet and hands of adults.

Size M: Best for applications to bigger muscle groups and medium sized joints, such as the arms, neck, wrist or elbows.

Size L: For the application to big muscle groups or joints such as the back, thigh, knee or ankle joint.

In our application tutorial, we will show you how to easily apply the Crosstapes by yourself. **To watch the video, simply scan the QR Code using the camera of your mobile phone:**


Our taping tips

Apply TRUETAPE® Crosstapes according to our instructions and follow these simple steps for the best taping results:

-  Thoroughly remove dirt, moisture, oils and lotions from the skin.
-  Trim excessive hair for better adhesion.
-  Apply tape min. 1 hour before starting an activity.
-  Apply Crosstapes pointing towards the head and feet. See more on next page.
-  To activate the adhesive, carefully rub the tape with your finger, starting from the inside moving towards the outside of the tape.
-  TRUETAPE® Crosstapes are waterproof. After showering simply pat the tape dry using a towel. Do not rub.
-  Don't rip it off! Gently remove tape from skin in the direction of hair growth.

Positioning of the tapes



For best wearing comfort and durability apply Crosstapes pointing approx. towards the head. See on left side below:

***Warning:**
The effect of Crosstapes has not been scientifically proven. Our statements are based on experiences and application reports. TRUETAPE® Crosstapes, along with the information and instructions we provide, are not a substitute for professional medical treatment and should not be used for self-diagnosis. Do not use Crosstapes in cases of acute inflammation, recent venous thrombosis, cancer, severe heart insufficiency, on body orifices, or on injured, irritated, or particularly sensitive skin. Crosstapes may cause skin irritation in rare cases. If you experience redness, irritation, pain, or itching, remove the Crosstapes immediately and carefully from the skin. During pregnancy, only use them in consultation with your physician. Do not use if you are allergic to any of the tape's components. Store tapes out of the reach of children.

Application guide: Pain points*

To apply TRUETAPE® Crosstapes to pain or trigger points, follow these steps:

-  Locate your pain point by palpating and moving the joints and muscles.
-  Select the optimal tape size. See also page 2.
-  Carefully remove the tape from the paper, touching only a small part of the adhesive.
-  Always apply the tape to stretched skin or bent joints.
-  Apply the tape with one tip pointing towards the head. See also page 4.
-  Rub Crosstape carefully with your finger from the inside to the outside.

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