Mamatest

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
Literatura

1. Batzer FR. Hormonal evaluation of early pregnancy. Fertil. Steril. 1980; 34(1): 1–13
2. Catti KJ, ML Dufau, JL Vaitukaitis Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyst. J. Clin. Endocrinol. Metab. 1975; 40(3): 537–540
3. Braunstein GD, J Rasor, H Danzer, D Adler, MC Wade Serum human chorionic gonadotropin levels throughout normal pregnancy. Am. J. Obstet. Gynecol. 1976; 126(6):678–681
4. Lenton EA, LM Neal, R Sulaiman Plasma concentration of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. Fertil. Steril. 1982;37(6): 773–778
5. Steier JA, P Bergsjo, OL Myking Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. Obstet. Gynecol.1984; 64(3): 391–394
6. Dawood MY, BB Savena, R Landesman Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma. Obstet. Gynecol. 1977; 50(2): 172–181
7. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross Ectopic production

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
Literatúra

1. Batzer FR. Hormonal evaluation of early pregnancy. Fertil. Steril. 1980; 34(1): 1–13
2. Catl KJ, ML Dufau, JL Vaitukaitis Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyst. J. Clin. Endocrinol. Metab. 1975; 40(3): 537–540
3. Braunstein GD, J Rasor, H Danzer, D Adler, ME Wade Serum human chorionic gonadotropin levels throughout normal pregnancy. Am. J. Obstet. Gynecol. 1976; 126(6): 678–681
4. Lenton EA, LM Neat, R Sulaiman Plasma concentration of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. Fertil. Steril. 1982; 37(6): 773–778
5. Steier JA, P Bergsjo, OT Myking Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. Obstet. Gynecol. 1984; 64(3): 391–394
6. Dawood MY, BB Savena, R Landesman Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma. Obstet. Gynecol. 1977; 50(2): 172–181
7. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross Ectopic production

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Hangzhou Biotest Biotech Co., Ltd., 17th Futai Road, Zhongtai Street, Yuhang District, 311121 Hangzhou, Zhejiang, China

Shanghai International Holding Corporation GmbH (Europe), Eiffelstraße 80, 20537 Hamburg, Germany

Distribútor: Green-Swan Pharmaceuticals CR, a.s., Tomičkova 2144/1, 148 00 Praha 4, Česko

www.gsklub.sk