ZuhauseTEST

DE

Gesunder Magen

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INSTRUCTIONS FOR USE

INTENDED USE

The **ZuhauseTEST Healthy Stomach** is a rapid chromatographic immunoassay for the qualitative detection of *H. pylori* antigens in human feces specimens, providing results in 10 minutes. The test utilizes antibodies specific for *H. pylori* antigens to selectively detect *H. pylori* antigens in human feces specimens.

INTRODUCTION

H. pylori is a small, spiral-shaped bacterium that lives in the surface of the stomach and duodenum. It is implicated in the etiology of a variety of gastrointestinal diseases, including duodenal and gastric ulcer, non-ulcer dyspepsia as well as active and chronic gastritis.^{1,2}

Both invasive and non-invasive methods are used to diagnose *H. pylori* infection in patients with symptoms of gastrointestinal disease. Specimen-dependent and costly invasive diagnostic methods include gastric or duodenal biopsy followed by urease testing (presumptive), culture, and/or histologic staining.³ A very common approach to the diagnosis of *H. pylori* infection is the serological identification of specific antibodies in infected patients. The main limitation of serology test is the inability to distinguish current and past infections. Antibody may be present in the patient's serum long after eradication of the organisms.⁴ HpSA (*H. pylori* Stool Antigen) testing is gaining popularity for diagnosis of *H. pylori* infection and also for monitoring the efficacy of the treatment of *H. pylori* infection. Studies have found that more than 90% of patients with duodenal ulcer and 80% of patients with gastric ulcer are infected with *H. pylori*.⁵

The **ZuhauseTEST Healthy Stomach** is a qualitative, lateral flow immunoassay for the detection of *H. pylori* antigens in human feces specimens. In this test, the membrane is pre-coated with anti-*H. pylori* antibodies on the test line region of the test. During testing, the specimen reacts with the particle coated with anti-*H. pylori* antibodies. The mixture migrates upward on the membrane by capillary action to react with anti-*H. pylori* antibodies on the membrane and generate a colored line. The presence of this colored line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

TEST CONTENTS

- 1 test cassette in a sealed pouch
- 1 sample collection applicator with sample dilution buffer
- 1 stool collection paper
- 1 instruction pamphlet

Additionally required:

- 1 timer

STORAGE AND STABILITY

The test cassette can be stored at room temperature or refrigerated (2-30°C). The test cassette is stable through the expiration date printed on the sealed pouch. The test cassette must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

TEST PREPARATION

Bring the test cassette and sample dilution buffer at a room temperature (15°C to 30°C) before performing the test.

Please read all the information in this package insert before performing the test.

TEST PERFORMANCE

1. Wash your hands with soap and rinse with clear water. To collect fecal specimens: The stool specimen should be collected in the stool collection paper. Please use the stool collection paper, avoiding contamination of the specimen by taking precautions that the specimen or side of paper containing specimen does not come in contact with any contaminating objects including toilet cleaners.
2. Take the sample collection applicator and open it by screwing the blue cap carefully. Be careful not to spill solution from the tube.


3. Randomly stab the specimen collection applicator into the fecal specimen in at least **3 different sites**. Do not scoop the fecal specimen.
4. Screw on and tighten the cap onto the specimen collection tube, then shake the specimen collection tube vigorously to mix the specimen and the extraction buffer.


5. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the test is performed immediately after opening the foil pouch.

6. Remove the white protective cap at the end of the tube. Using a piece of tissuepaper, break the tip of the tube. Be careful not to spill or spatter the sample mixture from tube.

7. Invert the specimen collection tube and transfer **2 full drops of the extracted specimen** to the specimen well (S) of the test cassette. Avoid trapping air bubbles in the specimen well (S). Then start the timer.

8. Read results at **10 minutes**. Do not read results after 20 minutes.


EVALUATION OF TEST RESULTS

To read the test results simply determine whether a line is present or absent at the control (C) position. It does not matter how strong or weak a Control line is. Staining of the control line (C) serves as an internal progress check. It confirms sufficient sample application and correct test performance.

POSITIVE TEST RESULT

The test result is **positive** if a **colored line appears in the control field (C) and a colored line is visible in the test field (T)**. This result means that *H. pylori* antigen is present in the stool and you should consult a physician.


NOTE: The intensity of the color in the test line region (T) will vary depending on the concentration of *H. pylori* antigen present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive. The final diagnosis must be confirmed by a physician.

NEGATIVE TEST RESULT

The test result is **negative** if a **colored line appears in the control field (C) and no colored line is visible in the test field (T)**.

This result means that the presence of the *H. pylori* antigen in feces was not detectable. If there are any other gastrointestinal symptoms, please, visit a physician for further medical investigation.

INVALID TEST RESULT

If there is **no Control line (C) or only a Test line (T)** in the result window, the test did not run correctly and the results are not valid.

Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The **ZuhauseTEST Healthy Stomach** has been evaluated with specimens obtained from a population of symptomatic and asymptomatic individuals. The result shows that the sensitivity of the **ZuhauseTEST Healthy Stomach** is 97.6% and the specificity is 97.9% relative to other rapid test.

Method	Other Rapid Test		Total Results
	Positive	Negative	
Zuhause-TEST Healthy Stomach	83	2	85
	2	93	95
Total Results	85	95	180

Relative Sensitivity: 97,6% (95%CI:*91,8%-99,7%)
Relative Specificity: 97,9% (95%CI:*92,6%-99,7%)
Overall accuracy: 97,8% (95%CI:*94,4%-99,4%)
*Confidence Intervals

Intra-Assay

Within-run precision has been determined by using 15 replicates of four specimens: negative, low titer positive, middle titer positive and high titer positive specimens. The specimens were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 15 independent assays on the same four specimens: negative, low titer positive, middle titer positive and high titer positive specimens. Three different lots of the **ZuhauseTEST Healthy**

Stomach have been tested using these specimens. The specimens were correctly identified >99% of the time.

Cross-reactivity

Cross reactivity with following organisms has been studied at 1,0*10⁶ organisms/mL. The following organisms were found negative when tested with the **ZuhauseTEST Healthy Stomach**: *Acinetobacter calcoaceticus*; *Acinetobacter spp*; *Branhamella catarrhalis*; *Candida albicans*; *Chlamydia trachomatis*; *Enterococcus faecium*; *E.coli*; *Enterococcus faecalis*; *Gardnerella vaginalis*; Group A *Streptococcus*; Group B *Streptococcus*; Group C *Streptococcus*; *Hemophilus influenzae*; *Klebsiella pneumoniae*; *Neisseria gonorrhoea*; *Neisseria meningitidis*; *Proteus mirabilis*; *Proteus vulgaris*; *Pseudomonas aeruginosa*; Rotavirus; *Salmonella choleraesuis*; *Staphylococcus aureus*; Adenovirus

Interfering Substances

The following potentially Interfering Substances were added to HPG negative and positive specimens. Ascorbic acid: 20 mg/dL; Oxalic acid: 60 mg/dL; Bilirubin: 100 mg/dL; Uric acid: 60 mg/dL; Aspirin: 20 mg/dL; Urea: 2000 mg/dL; Glucose: 2000 mg/dL; Caffeine: 40 mg/dL; Albumin: 2000 mg/dL.

BIBLIOGRAPHY

1. Marshall, BJ, McGechie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-444.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med.(1990), 322: 909-916.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-296.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

WARNINGS AND IMPORTANT INFORMATION

- For self-testing *in vitro* diagnostic use only. The test is for *in vitro* diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use the after expiration date.
- Keep out of the reach of children.
- The used test should be discarded according to local regulations.
- Store in a dry place at 2-30°C (36-86°F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- If you use a container to collect the stool samples, it must be clean.
- Follow the indicated time strictly.
- Use the test only once. Do not dismantle and touch the test window of the test

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RÉFÉRENCES

- Marshall, BJ, McGechie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-444.
- Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med.(1990), 322: 909-916.
- Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-296.
- Cutler AF. Testing for Helicobacter pylori in clinical practice. Am J. Med. 1996; 100:355-415.
- Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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