
lionelo

User manual
Instrukcja obsługi
Bedienungsanleitung
Инструкция по эксплуатации
Manuale d'uso
Manuel de l'Utilisateur
Manual de usuario
Handleiding
Naudojimo instrukcija
Návod k obsluze

Használati utasítás
Manual utilizare
Bruksanvisning
Bruksanvisning
Brugsvejledning
Käyttöopas
Manual de instruções
Οδηγίες χρήσης
Kasutusjuhend
Lietošanas instrukcija

BabyAir Flow

Nasal aspirator
Aspirator do nosa
Nasensauger
Назальный aspirатор
Aspiratore per naso
Aspirateur nasal
Aspirador nasal
Neusaspirator
Nosies aspiratorius
Aspirátor do nosu
Aspirátor az orrhoz
Aspirator nazal
Aspirator för näsan
Neseaspirator
Aspirator til næsen
Nenäaspiraattori
Aspirador nasal
Αναρροφητής μύτης
Ninaaspirator
Deguna aspirators

1


2

3


4


5

6


7


8



A




B



C



3

1

2

Dear Customer!

In case of any questions or comments on the purchased product, please contact us:
help@lionelo.com

Manufacturer:

BrandLine Group Sp. z o. o.
A. Kręlewskiego 1, 61-248 Poznań, Poland

IMPORTANT INFORMATION

1. Never leave the device unattended while connected to a power socket. The device may only be charged using a SELV charger with an output rating of 5V DC 1A.
2. This device may be used by persons with reduced physical, sensory or mental abilities, as well as by persons without relevant knowledge and experience in the use of the device, provided that they are supervised or have been previously instructed as to the safe use of the device and the risks involved.
3. This product is not a toy. Do not allow children to play with the device. Keep the device out of the reach of children.
4. Children must not perform cleaning or maintenance of the device.
5. Before using the product for the first time, please read the user manual carefully and keep it for future reference.
6. The product is intended solely for home use and not for medical or commercial purposes.
7. The product is a hygienic device designed for suctioning secretions and mucus to cleanse the nasal cavity of children from birth to 12 years of age, intended for home use.
8. If the child begins to feel unwell while using the device, stop using the device immediately and seek medical advice.
9. Do not use the product for purposes other than intended. The manufacturer bears no liability for damage caused by inappropriate use of the product.
10. Limit the use of the product to only a few seconds at a time in uninterrupted mode.
11. Clean the aspirator tips and secretion tank after each use.
12. To prevent pressure imbalances in the eardrums caused by vacuum pressure, infants and young children should close their mouth or swallow saliva while using this product.
13. Insert the aspirator tip into the nasal cavity gently, without applying excessive force.
14. Do not use the device if it begins to malfunction, overheats excessively, shows

discoloration or bulging, emits unusual sounds or smells, or if any other abnormal phenomena occur. In such cases, contact an authorized service center.

15. Do not place other objects on the device while it is in use.
16. Do not use or leave the device in direct sunlight.

BATTERY WARNINGS

1. **WARNING:** This product contains a button battery. Swallowing a button battery poses a risk of severe internal chemical burns.
2. **WARNING:** Dispose of used batteries immediately. Keep new and used batteries out of the reach of children. If you suspect that batteries have been swallowed or inserted into any part of the body, seek medical advice immediately.
3. Do not charge batteries that are not designed to be rechargeable.
4. Do not mix different types of batteries or combine new and used ones.
5. Insert batteries into the compartment according to the correct polarity.
6. Always remove an exhausted battery.
7. Do not short-circuit the supply terminals.

Contraindications

- Nosebleed
- Acute head or neck injury or trauma
- Coagulopathy or haemorrhagic disorders

Before using the device, please read the user manual carefully. Incorrect use of this product can pose health hazards and lead to the risks/complications listed below:

- Mechanical injury
- Nosebleed or internal bleeding
- Mucosal bleeding
- Discomfort and pain

LIST OF ELEMENTS (FIG. A)

1. Aspirator
2. Sponge for aspirator tip x5
3. Cleaning brush
4. Aspirator tube
5. 3 types of silicone nasal tip x3
6. Secretion tank
7. USB-C cable
8. Tweezers

DEVICE DESCRIPTION (FIG. B)

1. Battery level indicator
2. Suction power level indicator
3. Device operation indicator light
4. Suction ON/OFF button
5. Bedside lamp button
6. Night light
7. Suction power level control
8. Aspirator on/off switch

DESCRIPTION OF THE SECRETION TANK (FIG. C)

1. Tank light switch
2. Handle with a button for illuminating the nasal tip
3. Tank
4. Tank lock
5. Locking button
6. Secretion tank

OPERATING INSTRUCTIONS

Charging

Insert the USB-C power cable into the aspirator (A, Fig. 1) and connect the other end to the recommended charger.

The battery level indicator (B, Fig.1) will show:

- 3 bars – battery fully charged
- 2 bars – battery half charged
- 1 bar – connect to charger
- A flashing icon indicates that the battery level has dropped below 15%.

First use (Fig. 2)

1. Remove the secretion tank, the aspirator tube, and the appropriate silicone tip from the storage compartment in the aspirator.
2. Place the silicone tip in the secretion tank.
3. Connect the aspirator tube to the secretion tank and the aspirator. To start the device, press and hold the switch (8, Fig. B). Adjust the suction strength using the + and - buttons, then press the suction power button (4, Fig. B).
4. The aspirator and the secretion tank are equipped with a light function (5, Fig. B and 1, Fig. C).

Disassembly for cleaning (Fig. 3)

1. Disconnect the aspirator tube from both the aspirator and the secretion tank.
2. Remove the silicone nasal tip from the secretion tank.
3. Turn the secretion tank upside down, gently rotate the tank lid lock by 80 degrees, and then pull it upwards.
4. Press the locking button (5, Fig. C) and pull the base of the secretion tank.
5. Slide the handle with the nasal tip illumination button in the secretion tank to detach it (2, Fig. C).

Assembly after cleaning (Fig. 4)

1. Attach the secretion tank lock (B, Fig. 2) and then slide it into the base (A, Fig. 2). Turn it clockwise to lock.
2. Hold the base of the secretion tank, press the locking button (5, Fig. C), and gently insert it into the tank.

Sponge replacement (Fig. 5)

1. To replace the sponge, remove the lock from the base of the secretion tank (see: disassembly).
2. Insert a new sponge into the base of the tank, then reattach the secretion tank lock by turning it counterclockwise.

Replacing the button battery (Fig. 6)

Disconnect the transparent secretion tank. Loosen the screw securing the battery cover on the inner side. Remove both depleted batteries, then insert two new LR43 button batteries into the compartment, paying close attention to their polarity.

TROUBLESHOOTING

Weak suction / improper operation

- Check whether the silicone tip and other components are properly assembled.
- Inspect the sponge at the base of the secretion tank.
- Make sure the device is charged.

CLEANING AND MAINTENANCE

- Gently remove the aspirator tip and unscrew the secretion tank.
- Wash the silicone aspirator tip and secretion tank in warm water with a mild detergent.

Important! To clean the secretion tank, remove the part containing the button batteries (2, Fig. C).

- If necessary, clean the main unit with a damp cloth. Never submerge the main unit in water or any other liquids.
- Before reassembling the device, check to make sure that all parts are clean and dry.
- The aspirator tips and the secretion tank can be sterilised by steam or UV radiation.

The photos are for illustrative purposes only, the actual appearance of the products may differ from that shown in the photos.

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EN: The product complies with the requirements of European Union directives. In accordance with Directive 2012/19/EU, this product is subject to separate collection. The product should not be disposed of with municipal waste as it may pose a risk to the environment and human health. Return the used product to the recycling centre for electrical and electronic equipment.

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EN: The symbol placed on the product indicates that the batteries or accumulators used in it should not be disposed of with household waste. Such waste must be properly sorted in accordance with applicable regulations. Used batteries and accumulators must not be discarded in unsorted waste, as they may pose a risk to the environment and human health. Users should utilize available selective collection points, which enable their return, recycling, and safe disposal. For detailed information on local battery and accumulator recycling regulations, please contact your municipal office, local waste management facility, or landfill.

Note: Lithium batteries may pose a fire or explosion hazard, especially if damaged, overheated, or improperly disposed of. The chemicals contained within them may also be harmful to health and the environment.

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