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## EN SUNSIMED KA MEDICAL DEVICE

### INDICATIONS

Medical device with very high sun protection which contributes to prevent cancer lesions (excluding melanoma) and pre-cancer lesions such as actinic keratosis and UV induced DNA lesions. Developed specifically for sensitive skin vulnerable to sun exposure.

### KEY INGREDIENTS

- TriAsorB™** : Our 1<sup>st</sup> organic sun filter that protects from UVB, UVA radiation and high energy visible blue light.
- Pre-tocophéryl** : A powerful antioxidant for a re-enforced cellular protection.
- Hyaluronic acid**: Powerful hydrating agent.

### PROPERTIES

Its light texture and pump bottle eases product application and re application especially on top of the head. 8h Hydration\*\*.

Dry touch. Fragrance-free. High tolerance. Formulated to minimize risks of allergic reactions. Non comedogenic.

\*\*Clinical study: Hydration kinetic, single application, 23 subjects.

### DIRECTIONS FOR USE AND RECOMMENDED DOSES\*\*\*

Apply evenly each morning to the face and to sun-exposed areas throughout the day. Reapply in case of prolonged exposure, heavy perspiration or contact with water. Avoid direct contact with the sun's rays by wearing protective clothing (hat, sunglasses, T-shirt, etc.).

The illustrations hereafter can be used as a guide for estimating the number of doses (= pumps) needed for the different parts of the body (face, neck, neckline, head, forearms, hand).

\*\*\*Warning: the level of protection can be reduced if you do not use the recommended amount of product.

### PRECAUTIONS FOR USE

Do not use if you are allergic to any of the ingredients. If an allergic reaction occurs, consult a doctor. Avoid all direct contact with eyes in order to prevent irritation. In case of contact with eyes, rinse immediately with clean water. If an eye irritation persists, consult a doctor. Keep away from the reach of children. Close product tightly after use. Do not apply to mucous membranes. Do not apply to damaged skin. In case of concomitant treatment, seek your doctor's advice. It is not recommended to apply another product to the area requiring protection; this can reduce the efficacy of the preventive treatment. If you experience any undesirable effects related to the use of this product, see a doctor. If you have reason to believe that a serious incident occurred linked to the use of this product, please report it to your national authority and to the manufacturer.

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