



**MASKA KN95 FFP2**

**1 sztuka**

**SKŁADANA MASKA  
OCHRONNA  
KN 95 FFP2**



**PL**

1. Rozłóż maskę.
2. Umieść wewnętrzną stronę maski pod brodą i pociągnij obiema rękami za gumki.
3. Zabezpiecz gumki za uszami i ureguluj, aż będą wygodne.
4. Obiema rękami dopasuj nosek do nasady nosa, aż maska dobrze przylgnie. Sprawdź dopasowanie maski.

**ENG**

1. Open the mask.
2. Place the open side of the mask under the chin and pull on the elastic bands with both hands so that the nose-piece fits your face.
3. Secure the elastic bands behind your ears and adjust till comfortable.
4. With both hands, adjust the nose-piece to the bridge of your nose until the mask fits securely. Check the fit of the mask.



**Wytwórca**

**Zhejiang Benniao Technology Co., Ltd**

No. 159, Tangziyan Road, Changhe Street, Binjiang District,  
Hangzhou City, Zhejiang Province, China

**EC REP**

**Shanghai International Holding Corp. GmbH (Europe)**  
Eiffelstrasse 80, 20537 Hamburg, Niemcy

**Importer**

**PMT Sp. z o.o.**

ul. Tetmajera 79, 31-352 Kraków

[www.pmtmedical.pl](http://www.pmtmedical.pl)



REV 05/2020



5 907467 580487 >